

RESEARCH

Open Access



Spectrum subtraction as a complementary method for six resolution techniques resolving overlapping spectra; application to multicomponent veterinary formulation with greenness and whiteness assessment

Mahmoud G. Hagag¹ , Ahmed M. Hemdan¹ , Ahmed H. Nadim^{2*} , Samah S. Abbas² and Nesma M. Fahmy¹

Abstract

Mathematical filtration is an efficient tool to resolve the overlapping spectra of binary mixtures in zero or first order form. Herein, a comparative study was conducted between six economic, accurate and precise spectrophotometric methods for determination of Triclabendazole (TCB) and Levamisole HCl (LVM). Each component was resolved with minimum mathematical steps in its zero-order absorption spectrum by ratio subtraction, constant multiplication, and the recent factorized response method; coupled with spectrum subtraction. In addition, the mixture was resolved in its first derivative form by derivative subtraction, D^1 constant multiplication, and the recent D^1 factorized response method; coupled with spectrum subtraction. Results obtained were also compared to those obtained from constant value, concentration value, and derivative ratio methods. The linearity range was found to be either 1.0–10.0 $\mu\text{g/mL}$ or 2.0–20.0 $\mu\text{g/mL}$ for TCB, and 2.0–14.0 $\mu\text{g/mL}$ for LVM with LOD of 0.08 $\mu\text{g/mL}$ and 0.19 $\mu\text{g/mL}$, respectively. Validation of the proposed methods was performed according to VICH guidelines. Results obtained from the statistical data showed no significant difference regarding accuracy and precision compared to the reported methods. The developed spectrophotometric methods followed the principles of green analytical chemistry, in which the green assessment was done through four tools, called, National Environmental Methods Index (NEMI), Analytical Eco-Scale (AES), Green Analytical Procedure Index (GAPI) and Analytical greenness metric (AGREE). Also, a white assessment was performed using RGB model. The proposed methods could offer an economic alternative for the routine analysis of bulk materials and combined veterinary dosage form.

Keywords Green analytical chemistry, Spectrum subtraction, Factorized response method, Concentration value method, Levamisole HCl, Triclabendazole

*Correspondence:

Ahmed H. Nadim


ahmed.nagib@pharma.cu.edu.eg

Full list of author information is available at the end of the article



© The Author(s) 2023. **Open Access** This article is licensed under a Creative Commons Attribution 4.0 International License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit <http://creativecommons.org/licenses/by/4.0/>. The Creative Commons Public Domain Dedication waiver (<http://creativecommons.org/publicdomain/zero/1.0/>) applies to the data made available in this article, unless otherwise stated in a credit line to the data.

Graphical Abstract



Introduction

Triclabendazole (TCB), 5-Chloro-6-(2,3-dichlorophenoxy)-2-(methylsulfonyl)-1H-benzimidazole is utilized regularly beginning around 1983 to treat fascioliasis in the veterinary medicine causing inhibition of microtubule formation [1, 2]. It is considered as a powerful inhibitor of the synthesis of proteins used in goat, sheep, and cattle for the treatment of fascioliasis [3]. Although US Food and Drug Administration (FDA) has not supported the use of TCB in humans, however it's enlisted throughout many nations where fascioliasis is endemic [1, 4]. Its official Ph. Eur. method is potentiometric titration [2]. According to literature, it was reviewed in biological fluids, its metabolites presence, pharmaceutical dosage forms and pure form by spectrophotometric methods [5, 6] and chemometric methods [7, 8]. Literature also reveals its determination via HPLC [9–15] and LC–MS/MS [16–18]. Levamisole HCl (LVM), 2,3,5,6-Tetrahydro-6-phenylimidazo[2,1-b]thiazole monohydrochloride is utilized to treat parasitic worm infection. The drug host defenses by modulating cell-mediated immune responses. It restores depressed T-cell functions than stimulating response to above-normal levels. LVM is an anthelmintic with a broad spectrum of activity which is efficient against nematodes in both mature and immature stages, as well as

anthelmintic-resistant strains [19]. Its official method in Ph. Eur., USP and BP is a potentiometric titration method [2, 20, 21]. According to literature, it was reviewed in biological fluids, pharmaceutical dosage forms and pure form by spectrophotometric methods [22–24] and chemometric methods [8]. Literature also reveals its determination via HPLC [25–31], LC–MS/MS [32, 33] and GC [34]. Both drugs chemical structures are represented in Fig. 1.

TCB and LVM mixture has been utilized especially for the treatment of acute and subacute fasciolosis outbreaks as well as the essential handling of fascicle infections. There are only two reported spectrophotometric methods for the binary mixture analysis of TCB and LVM either using chemometric tools where sophisticated software was required or derivative method in which LOD and linearity could be significantly enhanced [7, 8]. Different spectrophotometric resolution methods can be applied. Ratio subtraction method coupled spectrum subtraction (RS–SS) and derivative subtraction method coupled spectrum subtraction (DS–SS) was used to get the less extended component X through the division of the mixture spectrum with a divisor of the more extended component Y' [35]. To get the more extended component (Y), spectrum subtraction (SS) is applied by subtraction of the




Fig. 1 Chemical structures of (a) Triclabendazole, b Levamisole HCl

X spectrum from the total mixture spectrum [36]. The same steps were repeated in DS-SS [36, 37] but on the first derivative spectrum mixture and utilizing a divisor spectrum in its first derivative, resulting in the first derivative spectra of both components allowing measurements to be done at its P_{\max} . Constant multiplication (CM-SS and D^1 CM-SS) method depended on determining the concentration of Y through the division of the mixture spectrum with a divisor of Y' (the extended component). At the extended region, the constant is measured and then multiplied by the mixture spectrum for getting the Y component (the more extended one) [37–39]. Also, spectrum subtraction method (SS) can be applied, where X could be determined by the subtraction of the Y spectrum from the total mixture spectrum [36]. The same steps were repeated in D^1 CM-SS [36, 37] but on the first derivative spectrum mixture and utilizing a divisor spectrum in its first derivative, resulting in the first derivative spectra of both components allowing measurements to be done at its P_{\max} . For factorized response spectrum (FC-SS), one response value (absorbance or derivative spectrum amplitude) which corresponds to component Y where X had no contribution ($\lambda_{\text{zero point}}$) using either its zero or derivative form has to be defined. The factorized response spectrum was measured by the division of (zero or derivative) every pure Y concentration by its recorded one at the selected wavelength. The spectrum of Y component obtained when the factorized spectrum is multiplied with this recorded response. By using spectrum subtraction (SS) method, X spectrum obtained by subtracting the Y spectrum from the total mixture spectrum [40–43]. For derivative ratio (DD^1), in which there were severely overlapping spectra of the D^0 absorption of the cited components. Then by using a divisor of X, the amplitude of

derivative ratio spectra mixture was directly proportional to the concentration of Y showing no X interference. Likewise, X can be determined by using a divisor of Y [44]. Constant value (CV) is one of the constant center method's two complementary steps [36, 38, 45] and it was utilized recently to analyze binary mixtures [39, 46]. It depended on the Y analysis in their binary mixture where X is less extended than Y. Through the division with a divisor spectrum, a constant from the plateau region parallel to the x-axis was obtained. The substitution process was done by utilizing its regression equation showing correlation between the amplitudes of ratio spectra Y/Y' at extended region. The corresponding concentration of Y at zero contribution of X, the Y concentration was determined. Concentration value (CNV) is a recent approach dealing with graphical manipulation of the spectra as the drug concentration value was determined directly indicating the actual concentration without any need of the substitution in its regression equation. It was performed through the division of D^0 mixture spectrum with the D^0 normalized spectrum (the more extended component). Then, the more extended component concentration was obtained from the constant determined from the plateau region parallel to the x-axis. [47, 48].

The aim of the current study was to establish a comparable study to develop six simple, economical, sensitive spectrophotometric methods which have the ability to resolve the spectral overlap of TCB and LVM in its dosage form. This work has been enhanced regarding the safety of humans and the environment due to the green ecofriendly methods which are safe for the planet also the sustainability ensured by white chemistry so it will not harm the future generations. This should offer a cost-effective alternative for the analysis of veterinary formulations.

Experimental

Device and software

A double-beam UV/Visible spectrophotometer model V-760 (Jasco, Japan) was used in spectrophotometric measurements. Data acquisition was done using Spectra manager[®] software JASCO corporation version 2 (Japan). The quartz cells used were obtained from Chromtech (UK) with 45 mm×12.5 mm×12.5 mm, 3.5 mL (1.0 cm)—2 optical faces. The absorption spectra of the blank and sample solutions were tested over 200–400 nm range.

Chemicals and reagents

Materials and reagents

Pure samples of TCB and LVM were provided from Pharma Swede company, Egypt with purity of 99.88 ± 0.84 and 100.32 ± 1.44 respectively, in accordance with Ph. Eur.

official methods [2]. Eradex Forte Suspension[®] each 1 mL of consists of 120 mg TCB and 75 mg of LVM was also obtained from Pharma Swede company, Egypt. Methanol and HCl of analytical grade were purchased from Piochem company, Egypt.

Standard stock solutions

Standard stock solutions were made to contain 500.0 µg/mL of each, TCB and LVM in methanolic HCl solution. Working standard solutions were prepared to contain 50.0 µg/mL of each, TCB and LVM in methanolic HCl solution. Throughout the lab work, standard stock solutions were kept in refrigerator (2–8 °C) for up to 5 days in glass-stoppered volumetric flasks.

Procedure

Spectral characteristics

Through scanning 200–400 nm range against methanolic HCl as a blank, the D^0 absorption spectra of 10.0 µg/mL TCB and 10.0 µg/mL LVM were obtained. Then, by spectra manager[®] software, the derivatized spectra were obtained and kept on the computer.

Factorized spectra of TCB

The D^0 TCB spectrum was divided by the recorded absorbance at 304 nm and saved on the computer. The D^1 TCB spectrum was divided by the recorded amplitude at 310 nm and saved on the computer.

Linearity and calibration graphs:

Portions equivalent to 10.0–200.0 µg TCB and 20.0–140.0 LVM were transferred from the working solutions (50.0 µg/mL) for TCB and LVM, into two volumetric flasks series of (10-mL) and completed to volume with methanolic HCl for TCB and LVM. The prepared standards spectra were scanned in the range of 200–400 nm against a blank then stored on the computer.

Calibration graphs for the zero order absorption spectra (D^0) of TCB and LVM Regression equations showing the linearity of relationships between the absorbance at λ_{\max} of the scanned spectra of TCB at 304 nm and 220 nm, and of LVM at 214 nm, versus the corresponding TCB or LVM concentrations were calculated.

Calibration graphs for the first derivative (D^1) spectra of TCB and LVM The first derivative spectra were determined using $\Delta\lambda=4$, and a scaling factor = 10. For obtaining the regression equations, calibration graphs were performed relating the peak amplitude (P-Zero) of the D^1 spectra of TCB at 310 nm and at $P_{\max-\min}$ ($P_{228-216}$) and

of LVM at P_{\max} ($P_{221\text{ nm}}$) versus the corresponding concentrations.

Calibration graphs for the first derivative of ratio spectra (DD^1) of TCB and LVM Regression equation relating amplitudes of the ratio spectra D^1 of TCB at $P_{\max-\min}$ ($P_{230-220\text{ nm}}$) and LVM P_{\max} ($P_{219\text{ nm}}$) versus the concentrations were calculated.


Calibration graphs for constant value method (CV) D^0 TCB spectra was divided by a divisor of normalized TCB. For computing the regression equation, the calibration graph was performed by plotting amplitude at the plateau region 290–306 nm versus the corresponding concentration of TCB.

Applying the spectrophotometric methods for determining TCB and LVM in laboratory prepared mixtures

TCB and LVM aliquots were transferred from their working solutions separately. Methanolic HCl was added to achieve the needed volume and obtain mixtures of various ratios of drugs under study and proceeded to determine the binary mixture concentration of TCB and LVM. Through the substitution of the absorbance of TCB at 304 nm in its corresponding regression equation, TCB concentration was obtained. Then, steps in scheme 1 were applied. Also, from the first derivative form of the lab mixture; for having TCB concentration, the peak amplitude at 310 nm was substituted in the corresponding regression equation. Then, the steps in scheme 2 were applied.

Ratio subtraction method coupled spectrum subtraction (RS-SS) and derivative subtraction method coupled spectrum subtraction (DS-SS) Through the division of the mixture D^0 absorption spectrum by TCB' divisor D^0 spectrum of 5.0 µg/mL, the constant (TCB/TCB) was obtained. This constant was subtracted from the D^0 binary mixture spectrum followed by multiplication with D^0 TCB' spectrum divisor of 5.0 µg/mL to get the less extended D^0 LVM spectrum to measure the absorbance as its λ_{\max} 214 nm. The D^0 TCB more extended spectrum was produced by subtraction of the resolved D^0 spectrum of LVM from the D^0 binary mixture spectrum, allowing measurement of the absorbance at its λ_{\max} 304 nm and 220 nm.

Through the division of the mixture D^1 derivative spectrum by TCB' divisor D^1 spectrum of 5.0 µg/mL, the constant (TCB/TCB) was obtained. This constant was subtracted from the D^1 binary mixture spectrum followed by the multiplication with D^1 TCB' spectrum divisor of 5.0 µg/mL to get the less extended D^1 LVM spectrum allowing measurement of the absorbance at



Scheme 1 Methods applied on D^0 for Triclabendazole and Levamisole HCl mixture


its P_{\max} (P 221 nm). The D^1 TCB more extended spectrum was produced by subtraction of the resolved D^1 spectrum of LVM from the D^1 binary mixture spectrum, allowing measurement of the amplitude at its P (310 nm) and $P_{\max-\min}$ (P 228–216).

Constant multiplication (CM-SS and D^1 CM-SS) Through the division of the mixture D^0 absorption spectrum by TCB' divisor D^0 spectrum of $5.0 \mu\text{g/mL}$, the constant (TCB/TCB') was obtained. The constant value in the extended region (290–306 nm) was reported and then multiplied with the $5.0 \mu\text{g/mL}$ TCB spectrum divisor to attain resolved D^0 more extended spectrum of TCB. Then, the absorbance at 220 nm in its corresponding regression equation was substituted. The D^0 less extended spectrum of LVM was obtained when the resolved D^0 TCB spectrum was subtracted from D^0 total mixture spectrum. Then, the absorbance at 214 nm in its corresponding regression equation was substituted.

Through the division of the mixture D^1 absorption spectrum by TCB' divisor D^1 spectrum of $5.0 \mu\text{g/}$

mL, the constant (TCB/TCB') was obtained. The constant value in the extended region (280–315 nm) was reported and then multiplied with the D^1 TCB spectrum divisor of $5.0 \mu\text{g/mL}$ to attain resolved D^1 more extended TCB spectrum. Then, the absorbance at $P_{\max-\min}$ (P 228–216) in its corresponding regression equation was substituted. From total mixture spectrum, when the resolved D^1 TCB spectrum was subtracted from the D^1 binary mixture spectrum, it led to the formation of the D^1 less extended spectrum of LVM. Then, the absorbance at P_{\max} (P 221) nm was substituted in its corresponding regression equation.

Factorized response spectrum (D^0 FC-SS and D^1 FC-SS) The D^0 TCB was determined through the multiplication of the peak amplitude at 304 nm with the TCB factorized spectrum λ_{\max} at 220 nm. Then subtracted from the mixture, to get LVM which was measured at 214 nm. The peak amplitude at 310 nm was multiplied



Scheme 2 Methods applied on D¹ for Triclabendazole and Levamisole HCl mixture

with the TCB factorized spectrum after derivatization to get its D¹ spectrum at P_{max-min} (P_{228-216 nm}). Then subtracted from the D¹ binary mixture spectrum to get the D¹ of LVM which have P_{max} (P_{221 nm}).

Derivative ratio (DD¹) of TCB and LVM The ratio spectra first derivative of TCB using LVM divisor of 4.0 µg/mL was determined, where the amplitude at P_{max-min} (P_{230-220 nm}) was substituted in the corresponding regres-

sion equation to obtain TCB concentration. The ratio spectra first derivative of LVM using TCB divisor of 3.0 µg/mL was determined and the amplitude at P_{max} (P_{219 nm}) was substituted in the corresponding regression equation to obtain LVM concentration.

Constant value (CV) The same constant derived from the plateau region 290–306 nm was substituted directly

in CV regression equation at 304 nm to obtain TCB concentration.

Concentration value (CNV) From D^0 spectrum of TCB, the peak amplitude was determined on the spectrum chart at plateau region 290–306 nm then R% was calculated.

Application to pharmaceutical formulation

One mL (1.045 is an average wt/mL) of Eradex Forte suspension[®] was transferred to volumetric flask (100-mL) and 50-mL of a solvent mixture of methanol: HCl (90:10, v/v) was added. The solution was then shaken for 5 min. The solution was sonicated for at least 40 min. and then filtered 3 times. The filtrate was completed with the same solvent mixture of methanolic HCl to the mark to get a stock solution with claimed concentrations of 1200.0 and 750.0 $\mu\text{g/mL}$ TCB and LVM, respectively. One mL of the stock solution was transferred into 50-mL volumetric flask and the volume was completed to the mark with methanolic HCl to have a final concentration of working solution 24.0 $\mu\text{g/mL}$ and 15.0 $\mu\text{g/mL}$ for TCB and LVM, respectively. As mentioned under analysis of laboratory prepared mixtures, the assay was done.

Results and discussion

In the previous years, several resolution techniques were proposed to mathematically filter the spectrum of each component in binary mixture in either its D^0 or D^1 absorption spectrum. The resolution techniques such as ratio subtraction or derivative subtraction could resolve the less extended component spectrum, while methods such CM and extended ratio subtraction (EXRS) could resolve the spectrum of the more extended component [46]. Combining two techniques such as RS/EXRS [47], or RS/CM [40] could resolve separately each component in the binary mixture. But as each method was composed of many mathematical steps, spectrum subtraction proved to be the best complementary method to any resolution technique able to resolve one component separately from the binary mixture, with the least manipulation steps and hence highest accuracy, eliminating any error that might appear from the choice of divisors and measurements at critical points.

As shown in Fig. 2a, TCB could be determined at its λ_{max} 304 nm where LVM had no contribution. But by applying resolution techniques that resulted in TCB zero order absorption, spectrum filtered alone allowed measurements to be done at its higher sensitivity peak 220 nm and eliminating the interference from LVM. Similarly,




Fig. 2 a The zero order absorption spectra of Triclabendazole and Levamisole HCl b The first derivative absorption spectra of Triclabendazole and Levamisole HCl

Fig. 2b shows that TCB could be determined directly at its P max 310 nm. With applying resolution techniques that resulted in the first derivative form of TCB alone allowed measurements to be conducted at Pmax-min ($P_{228-216 \text{ nm}}$) with higher sensitivity and eliminating the interference from LVM.

In this study, SS was coupled with six techniques, namely RS, DS, D^0 CM, D^1 CM, D^0 FC, and D^1 FC resolution techniques. A comparative study was established among all the proposed resolution techniques and the results were compared to DD^1 , CV, and CNV methods. The absorption spectrum of the two drugs under investigation showed high overlapping as indicated in Fig. 2. The TCB spectrum was more extended than LVM in the wavelength region from 290–306 nm. TCB could be determined directly at a λ_{max} 304 nm. As presented in Fig. 2, TCB had 2 peaks at 220 nm with linearity range 1.0–10.0 $\mu\text{g/mL}$ and at 304 nm with linearity 2.0–20.0 $\mu\text{g/mL}$. Since the peak at 220 nm had higher sensitivity, several spectrophotometric methods were conducted to allow measurements at this λ_{max} . Six resolutions techniques were conducted to obtain both TCB and LVM separately in their zero-order spectrum or first derivative absorption spectra, allowing measurement of each

drug at its robust λ_{max} or P_{max} with highest accuracy and precision.

Ratio subtraction (RS-SS) and constant multiplication (CM-SS)

The first resolution technique succeeded in producing the D^0 of LVM which is the less extended component was ratio subtraction (RS). This allowed measurements to be conducted at its λ_{max} 214 nm. By the application of spectrum subtraction (SS), the more extended D^0 TCB spectrum was obtained when the less extended D^0 of LVM was subtracted from the D^0 of the mixture spectrum, letting measurements of TCB to be applied at its λ_{max} 304 nm and 220 nm. Thus, RS/SS succeeded to obtain the spectra of each component alone in its D^0 absorption spectrum.

Constant multiplication was the resolution technique that succeeded in obtaining the more extended component D^0 spectrum of TCB, letting detections to be conducted at its λ_{max} 220 nm. Then by coupling this method with SS, the less extended D^0 LVM spectrum was obtained, where the resolved D^0 TCB spectrum was subtracted from the D^0 total mixture spectrum, allowing its measurement at its λ_{max} 214 nm. Thus, CM/SS succeeded




Fig. 3 a The zero order absorption spectrum of Triclabendazole and Levamisole HCl mixture, 10 $\mu\text{g/mL}$ of each overlaid with 5.0 $\mu\text{g/mL}$ triclabendazole as a divisor b The ratio spectrum of the mixture of Triclabendazole and Levamisole HCl divided by 5.0 $\mu\text{g/mL}$ triclabendazole as a divisor c The resolved D^0 absorption spectrum of Triclabendazole d The resolved D^0 absorption spectrum of Levamisole HCl




Fig. 4 **a** The zero-order absorption spectrum of Triclabendazole and Levamisole HCl mixture, 10 $\mu\text{g}/\text{mL}$ **b** The resolved D^0 absorption spectrum of Triclabendazole **c** The first order absorption spectrum of Triclabendazole and Levamisole HCl mixture, 10 $\mu\text{g}/\text{mL}$ **d** The resolved D^1 absorption spectrum of Triclabendazole

to obtain the spectra of each component alone in its zero-order absorption spectrum D^0 as indicated in Fig. 3.

Factorized spectrum (FC-SS)

D^0 Factorized resolution technique succeeded in having the D^0 of the more extended component TCB at its λ_{max} 220 nm. Coupling this method with SS succeeded in obtaining the D^0 of the less extended component LVM letting measurements to be conducted at its λ_{max} 214 nm. Thus, D^0 FC/SS succeeded to obtain the spectra of each component alone in its D^0 absorption spectrum as shown in Fig. 4.

Both RS and CM depended upon using a divisor of the extended component, while D^0 factorized resolution technique depended upon using a single wavelength where LVM had no contribution, and without the use of any divisors. Although RS/SS, D^0 CM/SS, and D^0 FC/SS ended up in resolving each component in its zero-order absorption spectrum, but D^0 FC/SS proved to be the most time effective technique, with minimum manipulation steps and no divisor. Comparing RS/SS with D^0 CM/

SS where both deal with using a divisor of the extended component, where a constant was acquired from the plateau region parallel to the x-axis. In CM, the constant was multiplied by the divisor for obtaining D^0 spectrum of TCB and coupling it to SS the LVM spectrum was obtained. In RS, the constant was first subtracted from the mixture's ratio spectrum then, multiplied with the divisor spectrum to get the less extended component LVM and coupling it to SS the TCB spectrum was obtained. The fewer steps in CM/SS made it superior to RS/SS.

Derivative subtraction (DS-SS), constant multiplication (D^1 CM-SS) and factorized spectrum (D^1 FC-SS)

The three techniques applied on the zero order absorption spectra were repeated on the first derivative spectra mixture showing each component first derivative spectrum alone where LVM was measured at its P_{max} 221 nm and TCB was measured at its $P_{\text{max-min}}$ (228–216

nm). The first derivative spectra of both components showed a wider extended region (280–315 nm) where LVM had no contribution. Derivative subtraction (DS) resulted in the first derivative of the less extended component LVM and when coupled to SS, the extended component D^1 spectrum of TCB was obtained. D^1 CM succeeded in obtaining the D^1 spectrum of the extended component TCB and when coupled to SS, the D^1 spectrum of the less extended component LVM was obtained. Both DS/SS, and D^1 CM/SS generated the first derivative spectrum of each component depending on a D^1 5.0 $\mu\text{g}/\text{mL}$ TCB divisor. But, D^1 CM/SS showed fewer manipulation steps. The third resolution technique D^1 FC/SS was superior to the previous two techniques and succeeded to obtain the first derivative spectrum of each component with the least manipulation steps.

Constant value (CV) and concentration value (CNV) methods

The same constant that resulted from the division of the lab mixture spectrum by a divisor of normalized

TCB' was either substituted in the corresponding regression equation to have TCB concentration by CV method or was directly represented the TCB concentration depending upon the CNV method. Due to the severely overlapping spectrum of both components in the wavelength region (200–250 nm), DD^1 resolution technique methods in which [44] using LVM divisor of 4.0 $\mu\text{g}/\text{mL}$, the LVM contribution in the lab mixture is cancelled when the ratio spectrum is derivatized, and the peak amplitude at $P_{\text{max-min}}$ ($P_{230-220 \text{ nm}}$) is directly proportional to the TCB concentration. While using 3.0 $\mu\text{g}/\text{mL}$ TCB as a divisor, the TCB contribution in the lab mixture is cancelled when the ratio spectrum is derivatized and the peak amplitude at P_{max} ($P_{219 \text{ nm}}$) is directly proportional to the LVM concentration as shown in Fig. 5.

Application to pharmaceutical formulation

This study was performed on Eradex Forte suspension[®] dosage form which contains two components of TCB and LVM. According to literature, TCB was soluble in many methanolic solvents. Methanolic HCl was proven to be




Fig. 5 **a** The first derivative spectra of ratio spectra of Levamisole HCl (2.0–14.0 $\mu\text{g}/\text{mL}$) using 3.0 $\mu\text{g}/\text{mL}$ triclabendazole as a divisor **b** The first derivative of ratio spectra of Triclabendazole (2.0–20.0 $\mu\text{g}/\text{mL}$) using 4.0 $\mu\text{g}/\text{mL}$ Levamisole HCl as a divisor

Table 1 Assay parameters and method validation obtained by applying the proposed spectrophotometric methods

Drug Name	TCB					LVM				
	D ⁰ (304 nm)	D ⁰ (220 nm)	D ¹ (310 nm)	D ¹ (P _{228-216 nm})	DD ¹ (P _{230-220 nm})	CNV	CV	D ⁰ (214 nm)	D ¹ (P _{221 nm})	DD ¹ (P _{219 nm})
Range (µg/mL)	2.0–20.0	1.0–10.0	2.0–20.0	1.0–10.0	2.0–20.0	...	2.0–20.0	2.0–14.0	0.0037	0.0099
Slope	0.0533	0.1322	0.0049	0.0088	0.0306	...	1.0043	0.0834	–0.0002	–0.0004
Intercept	0.0011	–0.0198	0.0001	–0.0002	–0.0031	...	0.0177	–0.003	98.79±0.94	100.54±0.37
Accuracy (R%±SD)	100.78±0.61	101.05±0.88	100.84±0.46	100.002±0.80	98.76±0.59	101.01±0.44	100.41±0.37	99.35±0.76	0.9998	0.9998
Correlation coefficient (r)	0.9998	1.0000	0.9999	0.9997	0.9999	...	0.9998	0.9997	0.365	0.203
Repeatability (RSD) ^a	0.176	0.522	0.334	0.410	0.585	0.100	0.184	0.537	0.383	0.769
Intermediate precision (RSD) ^b	0.554	1.159	0.386	0.531	0.821	0.216	0.259	0.565	0.19	0.23
LOD	0.28	0.08	0.26	0.14	0.52	...	0.29	0.27	0.58	0.70
LOQ	0.86	0.23	0.78	0.41	1.58	...	0.89	0.80		

RSD^a & RSD^b: relative standard deviation for the intra-day & inter-day precision, respectively (n = 3) using concentrations 4, 6, 9 µg/mL for TCB and 6, 9, 13 µg/mL for LVM

the best solvent which gave the most sensitive peak [5]. The literature had ensured that LVM was freely soluble in 0.1N HCl [22], so methanolic HCl was the solvent of choice. Many trials were performed to choose the best solvent ratio of methanol: HCl (50:50, 60:40, 70:30, 80:20, and 90:10). The ratio of methanol: HCl; 90:10 was chosen as the best ratio to extract the dosage form under study. Ultra-sonication time required for dosage form extraction was tried after 20, 25, 30, 35, and 40 min to get the best clear solution. The 40 min ultra-sonication proved to give satisfactory results. It should be noted that the proposed methods proved their sensitivity and selectivity in obtaining the spectra profile of TCB and LVM compared to literature [7, 8]. They spared the need of any sophisticated manipulation software. Also, the study had showed improved LOD for TCB and LVM of 0.08 and 0.19 $\mu\text{g}/\text{mL}$, respectively compared to 0.14 and 0.52 $\mu\text{g}/\text{mL}$ for the reported method. Additionally, these methods successfully proved their safety to human and environment. The methods are cost effective as well as can be utilized easily in the regular analysis of QC labs for drugs used in cattle breeding industry.

Validation

In compliance with VICH guidelines, methods validation was developed [49], where the parameters were presented in Table 1. They showed good accuracy, precision, and sensitivity. The specificity was evaluated through the binary mixtures analysis including a variety of drug ratios, and the methods were shown to be specific as indicated in Table 2. The developed methods were successfully used in the assay of cited drugs in Eradex Forte suspension[®]. Also, the results were statistically compared to reported methods as indicated in Table 3 and 4. Table 5 displayed the statistical results obtained by comparing the proposed methods to the reported method which show no significant difference.

Greenness assessment for the developed methods National environmental methods index (NEMI)

NEMI is regarded as one of the first methods for assessing green analytical chemistry. A circular diagram with four quadrants was constructed. Each one showed specific criteria including the usage of PBT (persistent, bioaccumulative and toxic) chemicals, consumption of hazardous chemicals, corrosiveness and the amount of waste generated [50, 51]. Three quadrants of the proposed methods showed green color as methanol doesn't exist either in the PBT list or in the TRI (Toxic Release Inventory) list, the pH of the study samples was non-corrosive (pH is not less than 2 and not higher than 12) and the amount of waste produced was less than 50 g/sample. While the only blank quadrant due to the presence of

methanol on the respective Resource Conservation and Recovery Acts lists as shown in Fig. 6.

Analytical eco-scale system (AES)

It was proposed by [52], depending on subtraction of the penalty points (PPs) which represent the harmful impacts of the method on the environment from 100. The AES assessment results showed that methanol have 6 PPs as it has 3 pictograms with the signal word "Danger" and uses less than 10 mL of sample volume for analysis (sample cuvette). Also, HCl which exists in very low concentration has only one pictogram with 2 PPs and its volume consumption during sample analysis (sample cuvette) is less than 10 mL. While for the waste in the instrumentation has 3 PPs because of 1–10 mL of waste with recycling by distillation. According to this assessment tool, when the score >75, it will be regarded as excellent green analysis. And when the final score is greater than 50, it will be regarded as acceptable green analysis. Inadequate green analysis is given for a final score less than 50. The AES assessment results are shown in Fig. 6 with excellent green analysis.

Green analytical procedure index (GAPI)

The GAPI assessment tool provides both semi-quantitative information and primary basic information. It has the ability to assess the whole analytical procedure [53]. The GAPI pictogram has 5 main sections with 15 detailed parameters. The color of the pictogram parts may be green, yellow, and red where the green color indicates a safe procedure while the red refers to non-ecofriendly operations. The developed spectrophotometric methods have 8 green shaded sections and 7 yellow shaded sections with no red shading. Figure 6 shows the pentagrams of the proposed method having low environmental influence.

Analytical greenness metric (AGREE)

Agree tool is a novel software which can be downloaded for green assessment [54]. It is considered a circle-colored pictogram with 12 principles of green analytical chemistry. Each section is colored from deep green (1) to deep red (0) depend on the environmental impact. The final score from 0 to 1 is displayed in the middle of that pictogram. All sections are colored deep green except Sect. 8 and 12 which are near deep green as the number of analytes in a single run are 2 with 20 samples per hour, methanol is considered highly inflammable respectively. While Sect. "Introduction", "Experimental" and "Conclusion" have a light green color due to the at-line analysis of the sampling procedure, 1 mL amount of sample and the low volume of waste respectively. Section 10 has yellow color because some of reagents used are bio-based.

Table 2 Determination of TCB and LVM in laboratory prepared mixtures and in Eradex Forte® suspension by the proposed methods

TCB:LVM ratio		TCB									
		D ⁰ (304 nm)		D ⁰ (220 nm)		D ¹ (310 nm)		D ¹ (P _{228-216 nm})		DD ¹ (P _{230-220 nm})	
R%±S.D	RS-SS	RS-SS	CM-SS	FC-SS	DS-SS	DS-SS	DS-SS	DS-SS	D ¹ CM-SS	D ¹ FC-SS	
	2	3	98.49±0.56	98.06±0.05	98.01±0.12	98.06±0.06	101.33±0.42	100.75±0.41	100.80±0.43	100.67±0.53	99.93±1.23
3	4	98.84±0.33	98.18±0.13	98.25±0.14	98.26±0.19	101.2±0.34	100.28±0.62	100.47±0.43	99.57±0.53	99.85±1.29	
1	1	99.86±1.94	98.29±1.99	100.58±1.78	101.09±1.94	101.18±0.98	100.43±0.69	100.99±0.59	100.81±0.49	99.99±1.07	
4.5	2.5	99.82±0.17	98.61±0.12	98.71±0.21	98.54±0.09	101.89±0.35	100.41±0.94	101.08±0.96	101.53±0.53	101.72±0.72	
*4	*2.5	98.57±0.75	101.06±1.97	100.3±1.76	101.12±1.94	100.31±1.53	100.49±0.62	101.81±0.66	99.34±0.53	99.53±0.99	
Eradex Forte® Suspension											
R%±S.D											
*4	*2.5	99.97±1.42	100.19±1.98	100.19±1.99	100.26±1.53	100.55±1.99	100.53±1.99	100.47±1.98	100.58±1.89	98.819±0.827	
TCB:LVM ratio		LVM									
		CV		CNV		D ⁰ (214 nm) FC-SS		D ¹ (P _{221 nm})		DD ¹ (P _{219 nm})	
R%±S.D	RS-SS	RS-SS	CM-SS	CM-SS	DS-SS	DS-SS	DS-SS	DS-SS	D ¹ CM-SS	D ¹ FC-SS	
	2	3	98.49±0.55	99.13±0.56	99.05±0.89	100.12±1.014	99.92±1.35	101.2±0.69	101.91±0.41	101.76±0.62	100.46±0.85
3	4	98.88±0.09	99.49±0.09	99.31±0.96	98.85±1.29	98.88±1.40	99.03±0.89	98.96±0.67	99.06±0.51	99.102±1.29	
1	1	98.93±1.83	99.61±1.83	101.7±0.50	101.2±0.59	101.46±0.97	100.77±0.77	100.51±0.97	100.39±1.39	100.14±1.44	
4.5	2.5	100.11±0.42	100.74±0.42	101.68±1.44	101.99±1.35	101.76±1.46	101.8±0.83	100.54±0.54	100.18±0.83	100.61±1.07	
*4	*2.5	100.47±1.19	101.13±1.19	101.93±1.07	101.48±1.26	101.44±0.96	100±1.95	100.72±1.65	100.3±1.65	100.21±0.94	
Eradex Forte® Suspension											
R%±S.D											
*4	*2.5	99.09±1.21	99.78±1.32	100.25±1.96	99.39±1.74	100.004±1.53	100.62±1.94	99.24±1.83	99.45±1.77	101.98±0.39	

*Dosage form ratio

Table 3 Statistical comparison between results obtained by the proposed methods and the reported methods for the determination of TCB and LVM in pure powder form

Parameters	TCB						LVM					
	Reported Method ^a [5]	D ⁰ (304 nm)	D ⁰ (220 nm)	D ¹ (310 nm)	D ¹ (P _{228-216 nm})	DD ¹ (P _{230-220 nm})	CV	CNV	Reported Method ^b [22]	D ⁰ (214 nm)	D ¹ (P _{221 nm})	DD ¹ (P _{219 nm})
Mean	99.88	99.92	99.73	100.39	99.81	99.85	99.88	100.48	100.32	100.17	99.35	100.10
S.D	0.84	0.90	0.87	0.83	0.85	1.49	1.03	1.35	1.44	1.26	0.82	1.15
n	6	6	6	6	6	6	6	6	6	6	6	6
Variance	0.71	0.81	0.76	0.69	0.72	2.22	1.06	1.82	2.07	1.59	0.67	1.32
F (5.05)	...	1.14 (5.05)	1.07 (5.05)	1.03 (5.05)	1.01 (5.05)	3.13 (5.05)	1.49 (5.05)	2.5896	...	1.30 (5.05)	3.09 (5.05)	1.57 (5.05)
Student's t ^c (2.228)	...	0.085 (2.228)	0.29 (2.228)	1.05 (2.228)	0.14 (2.228)	0.04 (2.228)	0.007 (2.228)	0.925 (2.228)	...	0.19 (2.228)	1.43 (2.228)	0.28 (2.228)

^a Spectrophotometric method using 0.1 M methanolic HCl as a solvent at wavelength of 305 nm

^b Spectrophotometric method using 0.1N HCl as a solvent at wavelength of 216 nm

^c The values are the corresponding tabulated values of t and F at p=0.05

Table 5 Results of ANOVA (single factor) for comparison of the proposed methods and the reported methods for determination of TCB and LVM in pure powder form


	Source of Variation	Sum of squares	Degree of freedom	Mean square	F	F crit
TCB^a [5]	Between groups	3.248	7	0.464	0.422	2.249
	Within groups	43.963	40	1.099		
	Total	47.211	47			
LVM^b [22]	Between groups	3.379	3	1.12	0.797	3.098
	Within groups	28.24	20	1.41		
	Total	31.62	23			
Comparison with the reported method for determining TCB and LVM in Eradex Forte [®] Suspension						
TCB^c [8]	Between groups	12.928	11	1.18	0.433	2.22
	Within groups	65.085	24	2.71		
	Total	78.013	35			
LVM^c [8]	Between groups	17.356	7	2.48	0.93	2.66
	Within groups	42.708	16	2.67		
	Total	60.06	23			

^a Spectrophotometric method using 0.1 M methanolic HCl as a solvent at wavelength of 305 nm

^b Spectrophotometric method using 0.1N HCl as a solvent at wavelength of 216 nm

^c Principal component regression (PCR) and partial least squares (PLS) chemometric methods

Green Assessment



White Assessment

Accuracy	101	Amount of reagents	66.6	Cost-effectiveness	100
Precision	91	Energy Intake	99	Time-effectiveness	98
Sensitivity	66.6	Other Hazards	100	Number of Analytes per hour	100
CS red 89		CS green 85		CS blue 99	
MB= 90.7% (Final color White)					

Fig. 6 Greenness and White assessment of the proposed spectrophotometric methods

Sect. “Results and discussion” has orange color as the at-line analytical device position. The proposed spectrophotometric methods AGREE score is (0.78) as shown in Fig. 6.

White analytical chemistry tools (WAC)

The purpose of WAC is to examine the analytical techniques from a variety of angles, such as their effectiveness, their effects on environment, and their sustainability from an economical point of view [55]. The recently established Red–Green–Blue (RGB) model, a quantitative evaluation tool that enables sustainability by measuring the whiteness of the analytical method, makes available algorithms for global evaluations of analytical methods. The analytical data refer to red color includes accuracy, precision, linearity range and LOD. While the green color involves amount of reagents, amount of waste, toxicity of chemicals, energy consumption and additional occupational risks. The Blue color contains cost and time effectiveness, sample material consumption, instrument service frequency and risk of random instrumental malfunctions. A free excel sheet is used in the evaluation of WAC using RGB algorithm with a result out of 100. The total score of the whiteness for the proposed spectrophotometric methods is 90.7 as shown in Fig. 6.

Conclusion

The three methods RS/SS, D⁰ CM/SS, and D⁰ FC/SS succeeded in obtaining each component spectrum alone in its D⁰ absorption spectrum with good accuracy and precision. The less extended component LVM was obtained by RS, and the extended component TCB was obtained by either D⁰ CM or D⁰ FC resolution techniques. Each time the resolved spectrum obtained was subtracted from the spectrum of the lab mixture to produce the D⁰ absorption spectrum of the other component by SS, that proved to be an excellent complementary method for each technique.

Similarly, the three methods DS/SS, D¹ CM/SS, and D¹ FC/SS succeeded to generate the spectrum of each component alone in its first derivative form with good accuracy and precision. The less extended component LVM was obtained by DS, and the extended component TCB was obtained by either D¹ CM or D¹ FC resolution techniques. Each time the resolved spectrum acquired was subtracted from the first derivative lab mixture spectrum to derive the first derivative spectrum of the other component by SS, that proved to be an excellent complementary method for each technique that allowed discarding the whole spectrum of the interfering substance.

The proposed methods showed their sensitivity and selectivity with few mathematical manipulation steps. Each spectrum was filtered alone to its D⁰ spectrum which performed as spectral profile of each cited drug allowed establishing spectra typical to the pure components present in Eradex Forte suspension[®]. In comparison to the CV method, there were very few manipulation steps in the concentration value method, as the concentration was shown graphically without substitution in regression equations. The proposed methods were determined with acceptable accuracy and precision at their maxima. The evaluation of the established procedures for the studied mixture’s greenness and whiteness confirmed its safety, sustainability and improved time and cost effectiveness. The greenness of analytical procedures was assessed using the most efficient assessment tools GAPI and AGREE, which provided integrated data on the entire utilized technique. Compared to GAPI, AGREE make analysis completed in a short amount of time by the freely downloadable software. GAPI offers visual examination and an assessment of all the phases of an analytical process. ESA uses numbers to refer for the greenness so they can be utilized together as method evaluation tools. NEMI is the simplest and least efficient greenness tool and should be used in conjunction with other reliable assessment tools. As a result of varying detail levels for each element assessed, the tools are different from one another. As a result, the analyst must decide which metric will be most useful. The WAC is a simple tool to evaluate various analytical techniques including the effectiveness of the techniques, environmental impact, and economic considerations through a provided Excel work sheet. They could be strongly used in QC laboratories to regularly analyze the studied drugs in their pure bulk powders or dosage forms without the need for any prior separation steps.

Abbreviations

TCB	Triclabendazole
LVM	Levamisole
SS	Spectrum subtraction
RS-SS	Ratio subtraction method coupled spectrum subtraction
DS-SS	And derivative subtraction method coupled spectrum subtraction
CM	Constant multiplication
FC	Factorized response
DD ¹	Derivative ratio
CV	Constant value
CNV	Concentration value
NEMI	National Environmental Methods Index
AES	Analytical Eco-scale System
GAPI	Green Analytical Procedure Index
AGREE	Analytical greenness metric
WAC	White Analytical Chemistry tools

Author contributions

MGH: methodology; investigation; formal analysis; visualization; writing—original draft. AMH: supervision; validation; resources; writing—review & editing. AHN: methodology; investigation; data curation; writing—review & editing. SSA: conceptualization; project administration; supervision; manuscript revision. NMF: methodology; validation; data analysis; writing—review & editing. All authors read and approved the final manuscript.

Funding

Open access funding provided by The Science, Technology & Innovation Funding Authority (STDF) in cooperation with The Egyptian Knowledge Bank (EKB). This research did not receive any specific grant from funding agencies in the public, commercial, or not-for-profit sectors.

Availability of data and materials

The datasets used and/or analysed during the current study are available from the corresponding author on reasonable request.

Declarations**Ethics approval and consent to participate**

Not applicable.

Consent for publication

Not applicable.

Competing interests

The authors declare that they have no competing interests.

Author details

¹Analytical Chemistry Department, Faculty of Pharmacy, Ahram Canadian University, 6 October City, 4th Industrial Zone, Banks Court Street, Giza, Egypt.

²Analytical Chemistry Department, Faculty of Pharmacy, Cairo University, El-Kasr El-Aini Street, Cairo 11562, Egypt.

Received: 5 April 2023 Accepted: 1 August 2023

Published online: 15 August 2023

References

- Robinson MW, Hoey EM, Fairweather I, Dalton JP, McGonigle S, Trudgett A. Characterisation of a beta-tubulin gene from the liver fluke. *Fasciola hepatica* Int J Parasitol. 2001;31:1264–8.
- COUNCIL OF EUROPE. European pharmacopoeia. Strasbourg, Council of Europe. 2020.
- Stitt AW, Fairweather I, Mackender RO. The effect of triclabendazole ("Fasinex") on protein synthesis by the liver fluke. *Fasciola hepatica* Int J Parasitol. 1995;25:421–9.
- Sweetman SC. Martindale the complete drug reference. 36th ed. London UK: The Pharmaceutical Press; 2009. p. 157.
- Shrivastava A, Kumar S, Jain A. Spectrophotometric method for quantitative determination of triclabendazole in bulk and pharmaceutical. *Chron Young Scient*. 2011;2(2):90–90.
- Ramadan NK, Mohamed AO, Shawky SE, Salem MY. Spectrophotometric determination of triclabendazole by acid-dye complexation method in bulk and pharmaceutical formulation. *J Appl Pharm Sci*. 2012;2(1):128.
- Dinç E, Pektaş G, Baleanu D. Continuous wavelet transform and derivative spectrophotometry for the quantitative spectral resolution of a mixture containing levamisole and triclabendazole in veterinary tablets. *Rev Anal Chem*. 2009;28(2):79–92.
- Pektaş G, Dinç E, Baleanu D. Spectrophotometric simultaneous determination of levamisole and triclabendazole in tablets by principal component regression and partial least squares chemometric methods. *Revista de Chimie-Bucharest-Original Edition*. 2008;59(2):154.
- Takeba K, Fujinuma K, Sakamoto M, Miyazaki T, Oka H, Itoh Y, Nakazawa H. simultaneous determination of triclabendazole and its sulphoxide and sulphone metabolites in bovine milk by high-performance liquid chromatography. *J Chromatogr A*. 2000;882(1–2):99–107.
- Rashed NS, Zayed S, Abdelazeem A, Fouad F. Development and validation of a green HPLC method for the analysis of clorsulon, albendazole, triclabendazole and ivermectin using monolithic column: assessment of the greenness of the proposed method. *Microchem J*. 2020;157: 105069.
- Mottier L, Moreno L, Alvarez L, Virkel G, Lanusse C. Measurement of triclabendazole and its metabolites in liver flukes: method development and full validation. *J Pharm Biomed Anal*. 2004;35(5):991–9.
- Denžić Lugomer M, Pavliček D, Bilandžić N, Majnarić D. Determination of benzimidazole residues and their metabolites in raw milk using high performance liquid chromatography–diode array detection. *Mljekarstvo Časopis Za Unaprjeđenje Proizvodnje I Prerade Mlijeka*. 2017;67(3):231–8.
- Nischal K, Somshekar B, Abhilekha P, Sharadamma K, Radhakrishna P. A simple RP-HPLC method for estimation of triclabendazole and ivermectin in a pharmaceutical suspension dosage form. *J Current Pharma Res*. 2011;1(4):306.
- Shurbaji M, Al Rub MHA, Saket MM, Qaisi AM, Salim ML, Abu-Nameh ES. Development and validation of a new HPLC-UV method for the simultaneous determination of triclabendazole and ivermectin B1a in a pharmaceutical formulation. *J AOAC Int*. 2010;93(6):1868–73.
- Alvinerie M, Galtier P. Assay of triclabendazole and its main metabolites in plasma by high-performance liquid chromatography. *J Chromatogr B Biomed Sci Appl*. 1986;374:409–14.
- Kikuchi H, Sakai T, Okura T, Nemoto S, Akiyama H. Total determination of triclabendazole and its metabolites in bovine tissues using liquid chromatography–tandem mass spectrometry. *J Chromatogr B*. 2019;1109:54–9.
- Cai C, Zhang L, Xue F, Qiu M, Zheng W. Simultaneous determination of triclabendazole and its metabolites in bovine and goat tissues by liquid chromatography–tandem mass spectrometry. *J Chromatogr B*. 2010;878(30):3106–12.
- Cai C, Xue F, Wang Z, Xiao S, Zhang L. Simultaneous determination of triclabendazole and its metabolites in bovine and goat fat tissues by liquid chromatography–tandem mass spectrometry. *Food Anal Methods*. 2012;5(6):1260–6.
- Amery WKP, Bruynseels JPJM. Levamisole, the story and the lessons. *Int J Immunopharmacol*. 1992;14:481–6.
- Pharmacopoeia U. USP 37-NF 32. In: *US Pharmacopoeial Convention United States Pharmacopoeial convention, Rockville*. 2014.
- Pharmacopoeia BH. The British Pharmacopoeia Commission Secretariat of the Medicines and Healthcare Products Regulatory Agency (MHRA). TSO (The Stationery Office): Publisher; 2019.
- Syed SM, Farooqui Z, Gawale K, Hamde D. A validated UV Spectroscopic method for determination of Levamisole HCl. *Int J Pharm Res Appl*. 2020;5:97–101.
- El-Didamony AM. Spectrophotometric determination of benzydamine HCl, levamisole HCl and mebeverine HCl through ion-pair complex formation with methyl orange. *Spectrochim Acta Part A Mol Biomol Spectrosc*. 2008;69(3):770–5.
- Abdel-Aziz O, Hussien EM, El Kosasy AM, Ahmed N. Comparative study of univariate spectrophotometry and multivariate calibration for the determination of levamisole hydrochloride and closantel sodium in a binary mixture. *J AOAC Int*. 2016;99(4):941–7.
- Abdallaheem A, Shantier S, Abureid I, Gadkariem EA. Stability-indicating high-performance liquid chromatographic determination of levamisole hydrochloride in bulk and dosage forms. *J Chem*. 2020;2020:6137897.
- El-Kholy H, Kempainen BW. Liquid chromatographic method with ultraviolet absorbance detection for measurement of levamisole in chicken tissues, eggs and plasma. *J Chromatogr B*. 2003;796(2):371–7.
- Ravisankar P, Rao GD. Development and validation of RP-HPLC method for determination of levamisole in bulk and dosage form. *Asian J Pharm Clin Res*. 2013;6(3):169–73.
- García JJ, Diez MJ, Sierra M, Terán MT. Determination of levamisole by HPLC in plasma samples in the presence of heparin and pentobarbital. *J Liq Chromatogr*. 1990;13(4):743–9.
- Cholifah S, Farina Kartinasari W, Indrayanto G. Simultaneous HPLC determination of levamisole hydrochloride and anhydrous niclosamide in veterinary powders, and its validation. *J Liq Chromatogr Relat Technol*. 2007;31(2):281–91.

30. Vandamme TF, Demoustier M, Rollmann B. Quantitation of levamisole in plasma using high performance liquid chromatography. *Eur J Drug Metab Pharmacokinet.* 1995;20(2):145–9.
31. Seid KH, Berhanu T, Asres K, Ashenef A. Development and validation of RP-HPLC method for the simultaneous determination of levamisole HCl and Oxytocin and its application in the assay of veterinary bolus formulations. *Curr Pharm Anal.* 2021;17(5):679–89.
32. De Baere S, Cherlet M, Croubels S, Baert K, De Backer P. Liquid chromatographic determination of levamisole in animal plasma: ultra-violet versus tandem mass spectrometric detection. *Anal Chim Acta.* 2003;483(1–2):215–24.
33. Jedziniak P, Szprengier-Juskiewicz T, Olejnik M. Determination of benzimidazoles and levamisole residues in milk by liquid chromatography–mass spectrometry: screening method development and validation. *J Chromatogr A.* 2009;1216(46):8165–72.
34. Schenck FJ, Podhorniak LV, Wagner R. A highly sensitive gas chromatographic determination of levamisole in milk. *Food Addit Contam.* 1998;15(4):411–4.
35. El-Bardicy MG, Lotfy HM, El-Sayed MA, El-Tarras MF. Smart stability-indicating spectrophotometric methods for determination of binary mixtures without prior separation. *J AOAC Int.* 2008;91(2):299–310.
36. Samir A, Lotfy HM, Salem H, Abdelkawy M. Development and validation of simultaneous spectrophotometric and TLC-spectrodensitometric methods for determination of beclomethasone dipropionate and salbutamol in combined dosage form. *Spectrochim Acta Part A Mol Biomol Spectrosc.* 2014;128:127–36.
37. Lotfy HM, Tawakkol SM, Fahmy NM, Shehata MA. Successive spectrophotometric resolution as a novel technique for the analysis of ternary mixtures of pharmaceuticals. *Spectrochim Acta Part A Mol Biomol Spectrosc.* 2014;121:313–23.
38. Lotfy HM. Determination of simvastatin and ezetimibe in combined tablet dosage forms by constant center spectrophotometric method. *Int J Pharm Pharm Sci.* 2012;4(4):673–9.
39. Lotfy HM, Saleh SS, Hassan NY, Elgizawy SM. Univariate versus multivariate spectrophotometric methods for simultaneous determination of complex binary mixtures with overlapped spectra: a comparative study. *Anal Chem Lett.* 2013;3(2):70–84.
40. Lotfy HM, Mohamed D, Elshahed MS. Novel univariate spectrophotometric determination of the recently released solid dosage form comprising dapagliflozin and saxagliptin via factorized response spectra: assessment of the average content and dosage unit uniformity of tablets. *Spectrochim Acta Part A Mol Biomol Spectrosc.* 2019;222: 117120.
41. Lotfy HM, Saleh SS. Investigating advanced approaches based on iso-absorptivity coefficient in unresolved spectral signals of binary mixtures. *J Anal Methods Chem.* 2019;2019:7924821.
42. Lotfy HM, Omran YR. Novel absorptivity centering method utilizing normalized and factorized spectra for analysis of mixtures with overlapping spectra in different matrices using built-in spectrophotometer software. *Spectrochim Acta Part A Mol Biomol Spectrosc.* 2018;200:167–78.
43. Lotfy HM, Saleh SS, El-Maraghy CM. Advanced approaches for the treatment and amplification of weak spectral signals produced by critical concentrations in white multicomponent systems. *Spectrochim Acta Part A Mol Biomol Spectrosc.* 2020;224: 117339.
44. Nevado JB, Cabanillas CG, Salinas F. Spectrophotometric resolution of ternary mixtures of salicylaldehyde, 3-hydroxybenzaldehyde and 4-hydroxybenzaldehyde by the derivative ratio spectrum-zero crossing method. *Talanta.* 1992;39(5):547–53.
45. Lotfy HM, Hegazy MAM. Simultaneous determination of some cholesterol-lowering drugs in their binary mixture by novel spectrophotometric methods. *Spectrochim Acta Part A Mol Biomol Spectrosc.* 2013;113:107–14.
46. Lotfy HM, Hagazy MA-M. Comparative study of novel spectrophotometric methods manipulating ratio spectra: an application on pharmaceutical ternary mixture of omeprazole, tinidazole and clarithromycin. *Spectrochimica Acta Part A Mol Biomol Spectrosc.* 2012;96:259–70.
47. Lotfy HM, Fayed YM, Tawakkol SM, Fahmy NM, Shehata MAE-A. Spectrophotometric determination for the binary mixture of clotrimazole and dexamethasone in pharmaceutical dosage form. *Anal Chem Lett.* 2017;7(1):30–42.
48. Lotfy HM, Fayed YM, Tawakkol SM, Fahmy NM, Shehata MAE-A. Evaluation of graphical and statistical representation of analytical signals of spectrophotometric methods. *Spectrochim Acta Part A Mol Biomol Spectrosc.* 2017;184:61–70.
49. Validation Methodology Guideline on Validation of Analytical Procedures: Methodology VICH Topic GL2, Step 7, CVMP/VICH/591/98, European Medicines Agency, London, UK. 1998.
50. Keith LH, Gron LU, Young JLCr. *Green Anal Methodol.* 2007;107(6):2695–708.
51. Tobiszewski M. Metrics for green analytical chemistry. *Anal Meth.* 2016;8(15):2993–9.
52. Gałuszka A, Migaszewski ZM, Konieczka P, Namieśnik J. Analytical Eco-Scale for assessing the greenness of analytical procedures. *Trends Anal Chem.* 2012;37:61–72.
53. Płotka-Wasyłka JJT. A new tool for the evaluation of the analytical procedure: green analytical procedure index. *Talanta.* 2018;181:204–9.
54. Pena-Pereira F, Wojnowski W, Tobiszewski MJ. AGREE—Analytical GREENness metric approach and software. *Anal Chem.* 2020;92(14):10076–82.
55. Nowak PM, Kościelniak PJA. What color is your method? Adaptation of the RGB additive color model to analytical method evaluation. *Anal Chem.* 2019;91(16):10343–52.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

